EndoVac[™]

Apical Negative Pressure Irrigation System

Instructions for Use

Contents

| | |
|---------------------|-----|
| English..... | 7 |
| Français | 11 |
| Deutsch..... | 17 |
| Español..... | 22 |
| Italiano | 28 |
| Português..... | 33 |
| Dansk | 38 |
| Nederlands | 43 |
| Norsk..... | 48 |
| Suomi..... | 53 |
| Svenska | 58 |
| Bulgarian | 63 |
| Magyar | 68 |
| Polski..... | 73 |
| Hrvatski jezik..... | 78 |
| Srpski | 82 |
| Eesti..... | 87 |
| Lietuviškai | 92 |
| Română..... | 97 |
| Türkçe | 102 |
| Latviešu..... | 106 |

Figure 1

Figure 2


Figure 3

Figure 4

Figure 5

Figure 6

Figure 7

2.

Figure 8

Figure 9

Figure 10



Figure 11

3.

Figure 12

Figure 13

Figure 14


Figure 15

4.

Figure 16

Figure 17

5.

Figure 18

Figure 19

Figure 20

Figure 21

6.

Figure 22

Figure 23

7.

Figure 24

Figure 25

How Supplied

The EndoVac™ Complete Installation Kit (973-3003) includes all major components*:

- | | |
|---------------------------------|---|
| A. Handpiece (2) | H. 3 cc Syringes for EDTA (6) |
| B. Fingerpiece (2) | I. MDT Evacuation Tubing (Blue Connectors) (2) |
| C. Multi-Port Adapter (1) | J. Handpiece/Fingerpiece Evacuation Tubing (White Connectors) (2) |
| D. MacroCannulas 25mm (20) | K. Instructions for Use (not shown) |
| E. MicroCannulas** 25mm (20) | L. Quick Setup Guide (not shown) |
| F. 20 cc Syringes for NaOCl (6) | |
| G. Master Delivery Tips MDT(20) | |

*Most components above can also be re-ordered separately.

**MicroCannulas can be reordered in 21, 25, 31 mm and 50 mm lengths.

Indications for Use

The EndoVac™ system is intended for the delivery and evacuation of endodontic irrigation solutions during root canal procedures.

Contraindication

None Known

Description of Device

The EndoVac™ System is intended for the delivery and evacuation of endodontic irrigation solutions during root canal procedures. EndoVac™ Apical Negative Pressure Irrigation System is used in root canal therapy and is comprised of four main components. The Multi-Port Adapter is used to connect the Master Delivery Tip to the dental office Hi-Vac system. The Master Delivery Tip provides a constant flow of irrigant that lowers the risk of overflow due to continuous suction. Additional tips include the MacroCannula which is used to remove coarse debris after instrumentation, and the MicroCannula which features microscopic evacuation holes allowing the irrigant to flow to the apical termination.

Warnings & Precautions

- Prior to EndoVac™ irrigation, the clinician must take every precaution necessary to ensure the **integrity of the rubber dam seal** and any potential communication from the pulp chamber to the oral cavity are absolutely secure.
- **Always protect the patient's eyes** with safety glasses and clothing from sodium hypochlorite splatter or spill.
- **Never place the MDT's delivery tip closer than five (5) millimeters from the coronal opening** of any pulp canal.
- **Never aim the MDT directly into the orifice** of the canal as this may create a positive force of irrigant into the canal.
- In order for the MicroCannula to be placed to full working length (WL), it requires a minimum canal preparation size of **at least a size #35 / .04 taper**, to full Working Length (WL) or at least a size #45 when a non-tapered instrument, such as LightSpeed®LSX™, is used.
- Always use the EndoVac™ system in the sequence described herein. Skipping or deviating from the following steps can cause the EndoVac™ cannulae to clog.
- Pre-fill various syringes: 20 cc syringes with sodium hypochlorite, 3 cc syringes with EDTA, and optional syringes with irrigants of the clinician's choice. Make sure no air bubbles are trapped in the pre-filled syringes, as this will cause uncontrolled irrigant extrusion after releasing the plunger pressure.

Adverse Reactions

None Known

Preclinical Considerations

- The EndoVac's fluid mechanics depend on an intact clinical crown with an access opening measuring at least 6-8 millimeters from cavosurface angle to the pulp floor. If the clinical crown is compromised, create a temporary crown using a composite material.
- **Final apical shape must be at least a #35 with a 4% taper.** For non-tapered preparations, it is recommended that LightSpeedLSX be used as a finishing file to enlarge to at least a size #45. This will provide adequate space for the MicroCannula to reach full WL and help ensure a cleaner canal.

1. Multi-Port Adapter and MDT Set-up in the Operatory

Figure 1. The Multi-Port Adapter (black arrow) is placed in a Hi-Vac connection as it sits in the evacuation rack. Ensure suction control valve (white arrow) is fully open.

Figure 2. Install Fingerpiece/Handpiece tubing set (white connectors) into both holes (black arrows) at the top of the Multi-Port Adapter. Loop the clear tubing behind the evacuation rack (white arrows).

Figure 3. Install the Master DeliveryTip (MDT with blue connector) onto the 20 cc syringe that holds NaOCl.

Figure 4. Firmly insert one end of the MDT tubing set (blue connectors) to the blue connector of the MDT (black arrow).

Figure 5. Firmly insert the opposite end of the MDT tubing set (blue connectors) to the lower port of the Multi-Port Adapter (black arrow).

Figure 6. Hang the MDT on the screw (black arrow) attached to the front of the Multi-Port Adapter.

Figure 7. The assistant or dentist now has convenient access to the MDT at all times.

2. MDT Use During Instrumentation

Overview

The MDT (Figure 8) provides both a constant source of fresh sodium hypochlorite from its metal delivery tip (Figure 9 - white arrow) into the pulp chamber and the immediate removal of any excess irrigant via the plastic evacuation hood (black arrow) surrounding the delivery tip. This dual action provides a method of maintaining a pulp chamber brim full with fresh sodium hypochlorite. The delivery tip is placed just inside the access opening while the evacuation hood remains on the outside. A stream of irrigant is directed from the delivery tip at an axial wall and never should it be aimed directly into a pulp canal orifice.

Reminder: It is possible to create positive pressure into the pulp canal, thus creating the risk of a sodium hypochlorite accident — if the clinician fails to follow these points: (a) always discharge the irrigant from the MDT into a natural or temporary intact access opening measuring at least six to eight (6-8) millimeters from cavosurface angle to pulp floor; (b) always direct the irrigant stream at an axial wall approximately 45 degrees from the pulp canal's axial plane in molars, 60 degrees in premolars, and 90 degrees in anteriors; and (c) never place the MDT's delivery tip closer than five (5) millimeters from the coronal opening of any pulp canal. (Figure 9)

Use

Instrumentation generates significant debris that must be evacuated from the pulp chamber via evacuation and fresh irrigant replacement. Use the MDT to deliver approximately one (1) cc of sodium hypochlorite at each instrument change in order to evacuate the debris and refresh the irrigant in the pulp chamber.

Figure 10. To avoid spillage of excess irrigant adjust the MDT tip and Evacuation Hood by bending it into a suitable position for each tooth.

Gripping the Syringe/MDT

Gripping the syringe plunger with the palm of the hand (white arrow) rather than with the thumb will reduce hand fatigue. (Figure 11). The needle component of the MDT was intentionally designed with a small internal diameter to prevent spills via a controlled flow. Some doctors or assistants do not have enough hand strength to comfortably express the NaOCl through a 20 cc syringe. In this case it is recommended to use 10 cc syringes and change the syringe more often.

3. MacroCannula Clinical Set-up

Figure 12. Firmly insert the MacroCannula into the Handpiece.

Figure 13. Remove the white connector closest to the front of the Multi-Port Adapter (arrow).

Figure 14. Firmly insert it into the back of the Handpiece (arrow).

Figure 15. Use the MDT to deliver irrigant into the pulp chamber.

Note: Check that the MacroCannula is suctioning properly. If not, check the Hi-Vac system and that the connections were made properly.

4. MacroCannula Clinical Use

Figure 16. After instrumentation, significant debris may remain on the pulp canal walls and in fin areas. This debris could clog the holes of the MicroCannula later in the procedure. Accordingly, the MDT supplies a copious volume of irrigant (white arrow) as the MacroCannula is worked up and down rapidly starting at the orifice (black arrow) past the mid-root area and into the apical third.

Figure 17. This vigorous up and down action lasts 30 seconds per canal (white arrows). Failure to properly perform this procedure may result in unnecessary clogging of the MicroCannula. If after 30 seconds debris continues to come out of the canal, continue until solution appears clear.

5. MicroCannula Clinical Set-up

Figure 18. Select the appropriate length MicroCannula (21, 25, 31 or 50 mm). Leave the cap on the MicroCannula to prevent damaging it during insertion into the Fingerpiece.

Figure 19. During insertion, the base of the MicroCannula is held as shown. Use **extremely** firm pressure when the two are mated. This mating pressure insures that the MicroCannula will not separate during use, especially in a curved canal.

Figure 20. Remove the cap of the MicroCannula. Clinical Note: Use the provided rubber stopper or a Sharpie® marker to indicate Working Length. In this picture, a Sharpie marker (left) is used to indicate 18 mm and the marker (right) represents 20 mm. It is suggested to mark the longest and shortest working lengths and interpolate the middle distance.

Figure 21. After the MicroCannula has been inserted into the Fingerpiece and length marked, the Fingerpiece replaces the Handpiece on the clear tubing (white connector).

6. MicroEvacuation

Figure 22. The length of the Fingerpiece is designed to match that of a file handle, making insertion into the pulp canal easy. The Fingerpiece is used to guide the MicroCannula into canal curvatures and leave the access opening readily available for irrigant delivery from the MDT.

The final phase of EndoVac™ irrigation is completed using the MicroCannula to flush and evacuate two separate irrigants, NaOCl and EDTA. The total sequence and time of each protocol is broken down into three cycles ("MicroCycles").

- First MicroCycle: Dissolves and removes organic debris from the pulp canal walls using 5.25 – 6% sodium hypochlorite.
- Second MicroCycle: Disassociates and removes the smear layer thereby exposing the dentinal tubules by using 15 – 17% EDTA.
- Third MicroCycle: Dissolves and removes the contents of the now exposed dentinal tubules, by again using 5.25 – 6% sodium hypochlorite.

Purging Hydrolysis Bubbles

During the instrumentation phase, the use of NaOCl produces gas bubbles — the hydrolysis of organic material. These microscopic gas bubbles adhere to the canal walls and must be evacuated. This evacuation is realized by delivering irrigant in 10 second intervals, then stopping briefly to purge the entire canal system of NaOCl and the consequential bubbles.

Figure 23. Once the Micro is placed at full Working Length, the clinician may remove her/his fingers from the Fingerpiece and proceed with irrigant delivery via the MDT.

7. MicroCycle Technique — Single-Canal

MicroCycle 1: NaOCl, MDT, 20 cc Syringe, MicroCannula, 30 seconds active, passive 60 seconds.

1. Place MicroCannula to full WL and deliver NaOCl from MDT for 10 seconds. Stop delivery and watch for MicroCannula to suction (PURGE) all NaOCl from the canal.
2. Repeat previous step again.
3. Again deliver NaOCl for 10 seconds, but instead of purging quickly remove the MicroCannula while the MDT continues to deliver NaOCl. This is referred to as "charging" the canal.
4. Let the NaOCl work for at least 60 seconds (passive wait) before proceeding to the next step.

MicroCycle 2: EDTA, MDT, 3 cc Syringe, MicroCannula, 10 seconds, passive 60 seconds.

1. Place MicroCannula to full WL and deliver EDTA from MDT for 10 seconds, leaving the canal filled (charged) with EDTA for 60 seconds.

MicroCycle 3: NaOCl, MDT, 20 cc Syringe, MicroCannula, 30 seconds active, passive 60 seconds.

1. Place MicroCannula to full WL and deliver NaOCl from MDT for 10 seconds. Stop delivery and watch for MicroCannula to suction (purge) all NaOCl from the canal.
2. Repeat previous step for another 10 seconds.
3. Repeat again, but after delivery for 10 seconds remove the MicroCannula immediately before removing the MDT, leaving the canal filled (charged) with NaOCl.
4. Let the NaOCl work for at least 60 seconds before proceeding to the next step. While waiting, fitting of a gutta-percha point is optional.

Dry Canal: MicroCannula, Paper points.

1. Place MicroCannula to full WL for at least 3 seconds.
2. Confirm canal is dry or continue drying with paper points.

8. Troubleshooting

Unclogging the Cannulae

Correct use of MacroCannula to remove the gross debris will minimize MicroCannula clogging. Once the clinician gains proficiency in clinical use, clogging becomes a non-issue and happens infrequently. If either cannula should become clogged, try unclogging it by attaching the back end of either the Fingerpiece or Handpiece onto a syringe filled with water. Push the plunger and in most cases the hole(s) are immediately cleared, as noted in (Figures 24 and 25).

Retreatment

In the case of retreatment, all pulp canals must be fully prepared, cleared and dried of gutta percha solvent as much as possible before proceeding to both Macro and Micro evacuation. Still, clogging may occur as this is the very nature of retreatment. When this occurs, unclog as noted above.

In anticipation of excessive clogging during retreatment, adequate clearing of the canal before evacuation reduces the chances of this happening.

Failure to Initially Draw Irrigant

Failure to achieve Micro flow can occur if the MicroCannula is taken through the foramen. If any of the Micro's twelve (12) holes extend beyond the apical foramen, they will not draw irrigant. If the clinician experiences this, withdraw the MicroCannula 0.5 mm at a time until flow starts.

9. Cleaning & Sterilization

| | |
|------------------------------------|---|
| WARNINGS | <ul style="list-style-type: none">• All re-usable components must be cleaned and sterilized prior to each use.• Failure to properly clean the components could lead to inadequate sterilization.• Only use cleaning solutions tested for efficacy and compatibility with the device / equipment being used.• Always observe all applicable legal and hygiene regulations for practice and / or hospital.• Always wear protective gloves, glasses and a mask when handling contaminated instruments. |
| Limitations on reprocessing | <ul style="list-style-type: none">• Repeated processing has minimal effect on the Aluminum Handpiece, Finger-Piece and Multi-Port Adapter; end of life is normally determined by wear and damage due to use.• For the Master Delivery Tip (MDT), it is recommended that the tip be replaced after 10 cycles in the autoclave. Replace the MDT earlier at any signs of wear or damage.• Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use. |
| INSTRUCTIONS | |
| Point of use: | <ul style="list-style-type: none">• Following use, components should be reprocessed as soon as is reasonably practical.• Components that may be processed for re-use include the Aluminum Handpiece, Aluminum Finger-Piece, Aluminum Multi-Port Adapter and the Master Delivery Tip.• Micro and MacroCannulas and evacuation tubing must be removed from autoclavable components before reprocessing. These items are single use only and should be discarded to prevent cross-contamination between patients.• Remove excess soil by wiping components with a disposable wipe and/or running under tepid (potable) water, checking for signs of visible contamination. Dry thoroughly with a disposable wipe. |

| | |
|-----------------------------------|---|
| Manual Cleaning | <p>Cleaning</p> <ul style="list-style-type: none"> • Use hospital grade enzymatic detergent. Follow manufacturer instructions for cleaning detergent concentration. • Soak device in detergent for 15 minutes. • Remove soil with cleaning detergent by scrubbing all external surfaces with a soft bristle brush until all visible soil has been removed. • Use a small diameter brush or pipe cleaner to clean cannulation holes. • All exposed surfaces should be free of visible soil. <p>Rinsing</p> <ul style="list-style-type: none"> • Immediately after cleaning, rinse device under warm running water for 3 minutes. • Pay special attention to rinse all holes and lumens. • Air dry in clean area. Blow lumens with clean air using filtered air source. |
| Maintenance and Inspection | <ul style="list-style-type: none"> • Inspect to ensure all visible contamination has been removed. • Check for damage and wear: corrosion, discoloration, cracking. • Discard damaged or corroded components. If not visibly clean, repeat the process or dispose of components. |
| Packaging | <ul style="list-style-type: none"> • Place product in a single FDA cleared sterilization pouch. • Ensure pouch is large enough to contain the components without stressing the seals. |
| Sterilization | <p>Pre-vacuum Steam Sterilization USA Full cycle: Wrapped 132°C (270°F) for 4 minutes, with a dry time of 20 minutes Gravity Steam Displacement Sterilization Full Cycle: Wrapped 121°C (250°F) for 30 minutes, with a dry time of 15 minutes EU Full Cycle: Unwrapped 134°C (273°F) for 3 minutes, with a dry time of 10 minutes</p> |
| Storage | To maintain sterility, instruments should remain wrapped until ready for use. |
| Additional Information | When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded. |

The instructions provided have been validated by the manufacturer of the medical device as being capable of preparing a medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to fully achieve the desired result. This normally requires validation and routine monitoring of the process.

10. Warranty

SybronEndo warrants certain components of the EndoVac™ system against defects in material and workmanship for 12 months from the date of purchase, subject to proper usage. During that twelve-month period, SybronEndo will replace:

- Master Delivery Tips, if received clogged.
- MacroCannulas, if received in unusable condition.
- MicroCannulas, if received in unusable condition.

Defects caused by misuse, neglect, accident or abuse are not covered by this warranty. SybronEndo assumes no liability resulting from improper use, damage or breakage due to misuse of these components by the purchaser. SybronEndo assumes no liability for damage to the EndoVac™ system components, injuries to patients or users or other problems resulting from improper use of accessories or other materials not supplied by SybronEndo.

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Kerr Italia S.r.l.
Kerr Australia Pty. Ltd.