ENGLISH

Presentations:

1 syringe with 2 g and 5 tips.

Indication:

Tinted Max Seal is a light-cured dental sealant used on pits and fissures of primary and permanent posterior teeth and complex anatomy in preventing caries.

Composition:

Bis-GMA, UDMA, HDDMA, sodium fluoride, calcium fluoride, BHT, camphorquinone, co-initiators and colloidal silicon dioxide.

How to use:

1. Check for the presence of occlusal contacts in the regions of pits and fissures. **2.** Perform absolute or relative isolation, ensuring there is no fluid contamination of the enamel to be treated. **3.** Perform the acid conditioning with 37% phosphoric acid (Maquira Acid Gel), for 30 seconds. Wash with plenty of water and then dry the surface. **4.** The conditioning site will be opaque and whitish. In case of contamination by fluids repeat the acid conditioning procedure. **5.** Apply the sealant to the pits and fissures, preserving the sites of occlusal contact. **6.** Light cure the region for 20 seconds. Light curing with power above 400 mW/cm² is recommended. Less powerful equipment will need more time for curing. **7.** Recheck the occlusal contacts.

Precautions:

Periodically verify the integrity of the material. If the material is totally removed, proceed as with the first application. If part of the material is removed, check the adhesion with an explorer on the remaining material, clean the site, perform acid etching, and apply the sealant only to the uncovered area.

Precautions:

For exclusive use by dentists;
Close the syringe immediately after use;
Products such as eugenol or clove oil inhibit polymerization;
Do not use the product if the validity has expired;
Keep out of reach of children;
Use safety goggles when handling the photopolymerizer.

Transport and Storage:

Always maintain the packaging well closed straight after use. It is advisable to store the product in a dry place, protected from the light and at temperatures between 5 and 30 °C. Do not store near products with an eugenol base.

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